

# VivaSight-EB



## Key Benefits

- Compatible with VivaSight-SL for a complete airway management solution
- Compatible with Ambu® aScope™ 3 Slim single-use bronchoscopes

## VivaSight-EB

### Single Lung Ventilation made simple and safe

VivaSight-EB is a single-use endobronchial blocker designated for Single Lung Ventilation (SLV). For visual guidance and positioning the VivaSight-EB can be used in combination with the Ambu® aScope™ 3 Slim single-use bronchoscope or the VivaSight-SL single lumen tube with integrated camera, depending on clinical preference.

The VivaSight-EB, consists of a sterile, single-use, “steerable” balloon-tipped catheter that is visually guided to a selected bronchial segment to affect balloon blockage of the distal

airway. The angled distal tip of VivaSight-EB can be adjusted to facilitate navigation to the bronchi. When used in conjunction with the VivaSight-SL continuous monitoring throughout the procedure ensures that dislocation can be easily detected and corrected.

For Single Lung Ventilation the Ambu visualisation portfolio includes VivaSight-DL, a double lumen tube with integrated camera, VivaSight-SL, a single lumen tube with integrated camera, VivaSight-EB, as well as the aScope 3 portfolio of single-use bronchoscopes.

# Specifications

VivaSight-EB	Size	Length
EBB09700	9Fr	700 mm
<p>General: VivaSight-EB features a blue silicone cuff.</p> <p>Usage: Single-use only. Sterile (ETO) unless package has been opened or damaged.</p> <p>Pack: Supplied individually packaged.</p> <p>Storage: store in a dry, cool and dark place.</p>		

## Ordering specifications

VivaSight-EB	Order number	Quantity/box
VivaSight-EB	EBB09700	5
Ambu® aScope™ 3	Order number	Quantity/box
Ambu® aScope™ 3 Slim	402001000	5
VivaSight-SL	Order number	Quantity/box
VivaSight-SL 7.0 mm	TVT70100	10
VivaSight-SL 7.5 mm	TVT75105	10
VivaSight-SL 8.0 mm	TVT80110	10
Accessories	Order number	Quantity/box
VivaSight Max 7" LCD monitor	SA00014	1
VivaSight Direct RCA cable for external monitors	SA00013	1



VivaSight-EB used with Ambu® aScope™ 3 Slim



VivaSight-EB used with VivaSight-SL

Ambu A/S  
Baltorpbakken 13  
DK-2750 Ballerup  
Denmark

VivaSight is distributed by Ambu A/S