

King Vision® aBlade™

Video laryngoscope



Key Benefits

Affordable

- The King Vision® aBlade is an excellent solution to solve the issue of limited availability of advanced visualisation solutions in the OR.
- With King Vision® you can improve patient care by having high performance visualisation capabilities for every intubation, at an affordable price.

Portable

- King Vision® is light weight, self-contained, battery operated and water resistant. King Vision® can be used continuously for more than 90 minutes making it highly portable.
- The aBlades are individually packaged. Each blade is disposable, eliminating concerns of cross contamination. King Vision® offers the choice of two blade types depending on preferred tube delivery technique.

Durable

- King Vision® and aBlades have been designed to be your primary tool for all intubations.
- Video output capability is compatible with external monitor and recording devices.
- The robust, full-colour, non-glare display can resist repeated cleaning. The aBlade video adapter is reusable, designed with materials that can withstand repeated use.

King Vision® video laryngoscope

Visionary performance in the palm of your hand. When it comes to visualising the airway for routine or difficult intubations, you need a tool that you can depend on each and every time. King Vision® is designed for indirect visualisation of the vocal cords which minimise the need for tissue manipulation and improve the ability to see anterior airways. Such reliable performance is a necessity, not an option. King Vision® comes with the choice of two blade types that allow users to keep their preferred technique in terms of tube delivery, offering an unique flexibility.

Portable

The King Vision® display is light weight, self-contained, battery operated, and water resistant. The reusable display and video adapter come packaged in a protective case. aBlades are individually packaged making King Vision® very portable.

Affordable

The disposable blades allow for economical use of the King Vision® for most of your intubations. Combined with the high performance visualisation capabilities this offers great value while also delivering excellent patient care. Take the guess-work out of your intubations. See the vocal cords and confirm placement of the endotracheal tube.

Durable

The King Vision® is designed to be your primary tool for intubations. The high-quality display comes with a 1-year warranty. The robust, full-colour, smudge-proof, non-glare display can resist repeated cleaning and normal use wear and tear.



Ordering specifications

Item no.	Description	Quantity
KVLAB3	Standard aBlade	20/CS
KVLAB3C	Channeled aBlade	20/CS
KVLVA34	aBlade Video Adapter 34	1/CS
KVLABKIT3	aBlade Kit	1/CS
KVIS01VA34	1 Display and 1 Video Adapter 34	1/CS
KVIS01ADCASE	aBlade case	1/CS

Part	Material
Blade	Polycarbonate, TPE (Thermoplastic Polyurethane Elastomer)
Video adapter	Polycarbonate/ABS, TPE (Thermoplastic Polyurethane Elastomer)
Display	Polycarbonate/ABS anti-reflective coating on display window

Sales Company
Ambu A/S
Baltorpbakken 13
DK - 2750 Ballerup
Denmark
Tel. +45 7225 2000
Fax +45 7225 2053
www.ambu.com

Manufacturer
King Systems
15011 Herriman Boulevard
Noblesville, IN 46060
Tel. 317 776 6823
Tel. 800 642 5464
Fax 317 776 6827
www.Kingsystems.com

Distributed by
Ambu

CE US: Rx only