

King Vision® has been designed with the intention to make a revolutionary series of high performance portable video laryngoscopes. The King Vision® combines the convenience of a durable reusable, video display with a disposable blade at a price that's affordable for all of your intubations, difficult and routine.

Learn more at www.airwaylearning.com

PRODUCT SPECIFICATIONS

REUSABLE DISPLAY		DISPOSABLE BLADES	
Batteries	3 AAA (Alkaline recommended)	Anti-fog	Coating on distal lens
Battery Life	> 90 min. (Battery status indicator flashes red when batteries need to be changed.)	AP Height of Blade	13 mm (non-channeled); 18 mm (channeled)
Computerized Power Management System	Auto shut off, Auto white balancing	Blade Length	17 cm
Video Aspect Ratio	4:3	Blade Width	26 mm (standard); 29 mm (channeled); 16 mm at distal tip
Video Port	RCA connection to monitor with cable	Camera Chip	CMOS
Video Refresh Rate	30 frames per second	Camera Resolution	640 x 480 VGA
Video Resolution	320 x 240 (QVGA)	Disposable Material	Polycarbonate / ABS
Video Screen	OLED Display	ET Size	6.0 mm - 8.0 mm (Channeled Blade)
Video Screen Size	6.1 cm / 2.4" diagonal	Light Source	White LED
Material	Polycarbonate / ABS anti-reflective coating on display window		

PRODUCT CODES

CODE	DESCRIPTION
KVL03	Standard Blade, Size 3
KVL03C	Channeled Blade, Size 3
KVIS01	Reusable Display
KVLKIT3	Vision Display, Case, AAA Batteries, 3 Channeled and 1 Standard Size 3 Vision Blades
KVCABL	Custom Video-out Cable (2.75 m / 9')

BETTER PATIENT CARE

The King Vision® is ergonomically designed to provide minimal lifting of soft tissue and impact on teeth.



King Vision® video laryngoscope

Leading the way in visualisation

Portable. Affordable. Durable.

Sales Company

Ambu A/S
Baltorpbakken 13
DK - 2750 Ballerup
Denmark
Tel. +45 7225 2000
Fax +45 7225 2053
www.ambu.com

Manufacturer

King Systems
15011 Herriman Boulevard
Noblesville, IN 46060
Tel. 317 776 6823
Tel. 800 642 5464
Fax 317 776 6827
www.Kingsystems.com

Distributed by
Ambu



The growing challenge in healthcare around the world fuels our work to release the power of single-use devices and transform ideas into breakthrough solutions that work for life.

King Vision® video laryngoscope

Visionary performance in the palm of your hand

- Advanced imaging technology and unparalleled ergonomic design
- Choice of two blade types standard and channeled depending on the complexity of the airway
- Use the King Vision® for all your intubations



AFFORDABLE

- Disposable blades allow economical use of the King Vision® video laryngoscope for all intubations.
- With the King Vision® you can improve patient care by having high performance visualisation capabilities for every intubation, at an affordable price.



PORTABLE

- The King Vision® is light weight, self-contained, battery operated and water resistant.
- The King Vision® can be used continuously for more than 1 hour making it highly portable.
- Blades are individually packaged. Each blade is disposable, eliminating concerns on cross contamination.



DURABLE

- Designed to be your primary tool for all intubations.
- Video output capability is compatible with external monitors and recording devices.
- The robust, full-colour, non-glare display can resist repeated cleaning and normal use wear and tear.



REUSABLE
DISPLAY



DISPOSABLE
STANDARD
BLADE



DISPOSABLE
CHANNELED
BLADE