

# KingMask

anaesthesia breathing mask



## Key Benefits

### Market Leading Design

- Cushion is extremely pliable and tacky to the touch, for an excellent face seal
- Ribbed on top of the cushion, for a no-slip grip
- Anatomically shaped cuff enabling a tight seal with minimum applied pressure
- Flexible, self-expanding dome
- Clear dome for easy observation of the patient's condition
- Latex free
- DEHP free material

## KingMask

Because of the innovative quality design of the King Mask, it is the market leading anaesthesia mask design available today. The King Mask has an exceptionally thin and tacky cushion which improves the ability to create a ventilation seal on a variety of patients. The design with rounded edges secures a smooth grip and high user conformity.

The King Mask is also available in a scented version with Fresh Scent, Bubble Gum, Strawberry and Cherry scents. The Sweet Dreams™ masks release soothing aroma from natural oils in mask cushion. This can have a calming effect during induction, especially for children<sup>1</sup>.

### Five Scents

- Unscented
- Fresh Scent (mint)
- Bubble Gum
- Strawberry
- Cherry

### Seven Sizes

- Neonate
- Infant
- Toddler
- Child / Small Adult
- Medium Adult
- Large Adult
- X-Large Adult

<sup>1</sup>Carson KD et al. Quantitative assessment of patient preference for anaesthesia face mask. Irish Journal of Medical Science 1995;164 suppl. 14

# KingMask

anaesthesia breathing mask

Mask	Neonate (1)	Infant (2)	Toddler (3)	Child /Sm Adult (4)	Md Adult (5)	Lg Adult (6)	XL Adult (round)
Original (unscented)	1012 / 1015	1022 / 1025	1032 / 1035	1042 / 1045	1055	1065	1075
Fresh Scent (mint)	2012 / 2015	2022 / 2025	2032 / 2035	2042 / 2045	2052 / 2055	2062 / 2065	2072 / 2075
Bubble Gum		1122 / 1125	1132 / 1135	1142 / 1145			
Strawberry		1222 / 1225	1232 / 1235	1242 / 1245			
Cherry		1322 / 1325	1332 / 1335	1342 / 1345			

Part #'s ending in 2 = 20 per case • Part #'s ending in 5 = 50 per case



Unscented



Fresh Scent



Bubble Gum



Strawberry



Cherry

## Materials

Part	Material
Cushion	PVC Polyvinyl Chloride (DEHP free)
Crown	Polycarbonate
Hook ring	Polycarbonate
Valve	PVC Polyvinyl Chloride (DEHP free)
Valve Core	Polyester
Spring	Stainless steel

## Product weight

Size	Weight
1	16.6g
2	21.7g
3	31g
4	40g
5	52.6g
6	67.7g
XL Adult	72g

### Sales Company

**Ambu A/S**  
Baltorpbakken 13  
DK - 2750 Ballerup  
Denmark  
Tel. +45 7225 2000  
Fax +45 7225 2053  
www.ambu.com

### Manufacturer

**King Systems**  
15011 Herriman Boulevard  
Noblesville, IN 46060  
Tel. 317 776 6823  
Tel. 800 642 5464  
Fax 317 776 6827  
www.Kingsystems.com



US: Rx only