

At a higher level of practicality

The single use concept combined with its user-friendly lightweight design and new camera technology makes the aScope 2 a unique alternative to reusable scopes. After use it can just be discarded as with any other disposables used. In short, the Ambu® aScope™ 2 provides greater convenience for everyone involved.

Eliminating the risk of cross contamination

Flexible scopes are difficult to clean and disinfect because of their long and narrow lumens. Guidelines and recommendations should be strictly followed but are often not, as it is time consuming and costly.^{1,2,3} With the increasing risk of resistant bacterial strains, the need for new and smarter solutions such as the aScope 2 is evident.

With low investment, no cleaning and no repair

Reusable scopes have high initial investment costs and high maintenance costs as they are repeatedly damaged despite staff training. As a single use device, Ambu® aScope™ 2 is a cost-effective solution for departments and clinics of all sizes.

Main benefits of Ambu® aScope™ 2:

- Availability: Always there when you need it
- Eliminates the risk of cross contamination
- Cost effective: no cleaning and repair costs.
- Sterile and ready for use after connection to the Ambu® aScope™ Monitor.
- ClearLens design with optimised image quality.
- A powerful solution for monitoring during Percutaneous Dilational Tracheostomy.



aScope 2 for intubation

Try aScope 2 if you are tired of limited access to reusable scopes. Whether a difficult airway is anticipated or unanticipated aScope 2 is ready to use whenever you need, and it is affordable in even remote elective ORs

aScope 2 for PDT

The use of endoscopic guidance during PDT is increasing and is recommended for routine use⁴. The Ambu® aScope™ 2 offers a clear endoscopic view⁵ and allows the operator to see and verify the accurate placement of the guidewire and needle prior to tracheal dilatation as well as the placement of the tracheostomy tube.



1) Daniel J Weber, MD, MPH; William A. Rutala, PhD, MPH, Lessons from Outbreaks Associated With Bronchoscopy, Infection Control and Hospital Epidemiology, July 2001, Vol.22 No.7 pp403-408 2) Jennifer Schraag, Scope-related Infections, The same challenges have remained for decades, but have the processes improved? May 2007, www.infectioncontrolday.com 3) A. Norris, P. Wakeling, M. Wiles, R. McCahon, M. Bennett, Costs associated with fibrescope inventory: what role for a disposable scope?, DAS 2010 (Poster) 4) Kost K.M., Endoscopic percutaneous dilatational tracheotomy: a prospective evaluation of 500 consecutive cases. Laryngoscope. 2005 Oct;115(10Pt2):1-30 5) C Gernoth, H.V. Genzwuerker, Euroanaesthesia 2010, Helsinki, Finland – June 12-15, 2010

US: RX only

Ambu A/S
Baltorbakken 13
DK - 2750 Ballerup
Tlf.: +45 7225 2000
Mail: hello@ambu.com
www.ambu.com

Ambu
Ideas that work for life

The growing challenge in healthcare around the world fuels our work to release the power of single-use devices and transform ideas into breakthrough solutions that work for life.

Ambu® aScope™ 2
A Unique Single Use, Flexible Intubation Scope

Ambu® aScope™ 2

leading the way

The single use flexible intubation scope that guides difficult airway intubations and PDT procedures.

- Always available
- No cross contamination
- No cleaning or repair costs



Your Idea

"Why not a single use flexible intubation scope?"

When asking for the flexible scope you and your colleagues are often faced with answers like "The scope is being used somewhere else, it's being cleaned, it's submitted to maintenance, it's out of order, etc." This made you approach us with an apparently impossible thought: "Why can't we have a disposable scope ready to use straight from the box?"



Our work

Intubation and PDT made easy – on every occasion

Utilising advanced video technology and profiting from our networks vast expertise we carried out a vision: a flexible single use intubation scope that once and for all solves the problem of availability, simplifies set-up procedures and eliminates the need for cumbersome cleaning procedures – at a highly competitive price



Everyone's life

Enabling you to provide the best care

Difficult airways happen whether we plan them or not. With the aScope 2 you now have an accurate, affordable and convenient way to handle these critical situations with no risk of cross contamination. In helping you advance both patient care and safety, aScope 2 is a giant step in the right direction.

